



# IPS Empress® CAD RESTORATIONS

BY STRAUMANN® CARES®



COMMITTED TO  
**SIMPLY DOING MORE**  
FOR DENTAL PROFESSIONALS

## FEATURES AND BENEFITS

Straumann presents the IPS Empress® CAD restorations by Straumann® CARES®.

### **ALL-CERAMIC RESTORATIONS FOR OUTSTANDING NATURAL ESTHETICS**

- Material scatters light like natural enamel
- Multi-polychromatic material providing smooth color transition and fluorescence gradient
- Precise tooth design supported by Straumann® CARES® Visual software<sup>1</sup>, e.g. Ivoclar Vivadent AG tooth database

### **SIMPLICITY FOR EASY HANDLING**

- User-friendly and flexible software
- Staining or cut-back processing technique, depending on the application
- Adhesive cementation

### **WELL-ESTABLISHED PRODUCT DESIGNED FOR PREDICTABLE RESTORATIONS**

- 20 years of clinical experience with leucite glass-ceramics
- Straumann's industrial milling center utilizes cutting-edge technology designed to deliver outstanding quality

IPS Empress® CAD is a proven leucite-reinforced glass-ceramic from Ivoclar Vivadent AG, an established specialist in ceramic materials and final restorations. Straumann® CARES® CAD/CAM offers IPS Empress® CAD restorations of high quality due to Straumann's continuous investment in cutting-edge production technologies.

<sup>1</sup> Version 6.0 or higher.  
Some products may require regulatory approval and may not be available in all markets.

## APPLICATIONS AT A GLANCE

With the IPS Empress® CAD restorations, Straumann® CARES® CAD/CAM offers an outstanding solution for:

- Full contour inlays, onlays, veneers, partial crowns and crowns
- Cut-back veneers and crowns<sup>2</sup>

Depending on the processing technique and application, the following recommendation can be made to select the proper material:

Indications	Translucency level			Processing techniques	
	HT	LT	Multi	Staining	Cut-back
Inlays	●			●	
Onlays	●			●	
Veneers	●	●	●	●	●
Partial crowns		●	●	●	
Crowns		●	●	●	●

### Benefit from a leading material offering and application range with Straumann® CARES® CAD/CAM.

For further information visit our website or contact your local Straumann subsidiary.

<sup>2</sup> Applications may not be available in all markets.

IPS Empress® is a registered trademark of Ivoclar Vivadent AG, Liechtenstein.

© Institut Straumann AG, 2012. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates. All rights reserved.